

DR Liberty



*System installs in
under 1 hour*

Display

- ▶ 12.1" TN LED backlight digitizer wide screen (WXGA) display with dual touch

System Storage

- ▶ 128GB Solid State Drive (SSD)

DR Liberty Specifications

Wireless Flat Panel Specifications

- ▶ Wireless Cassette based Detector
- ▶ Image Size 35x43cm (14" x 17")
- ▶ Wireless Standard 802.11n
- ▶ 8.5 lbs (3.9kg)
- ▶ 139m Pixel pitch
- ▶ 70-90 images on a single charge
- ▶ Utilized with standard grids
- ▶ Liberty Software by INFIMED
- ▶ 24/7 Technical support team

Console Specifications

Processor / Cache

- ▶ Intel® Core™ i7-640UM vPro™ Processor 1.2 GHz with Intel® Smart Cache (4MB of L3 cache, 2 Cores/4 Threads)

Operating System

- ▶ Genuine Windows® 7 Professional (32-bit)

Durability Features

- ▶ Water, dust and splash resistant
- ▶ Shock-mounted display and hard drive
- ▶ Rubberized coated bottom housing for shock dampening

Batteries

- ▶ Lithium-ion battery with 30WHr capacity
- ▶ One battery capacity: ~3.3 hours
- ▶ Two battery capacity: ~7 hours

Warranty

- ▶ Standard 1-year field-ready warranty
- ▶ Multi-year options available



First Source, Inc.
1100 Iverleigh Trail
Charlotte, NC 28270
1-800-349-5980
www.fsmedicalimaging.com